



Product Datasheet

NEK9 (Phospho-T210) antibody (orb315630)


Description	Rabbit polyclonal antibody to NEK9 (Phospho-T210)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T210 of human NEK9 protein. The exact sequence is proprietary.
Target	NEK9
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8TD19, Q8K1R7
Entrez	217718, 91754
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of NEK9 (Phospho-T...



Immunohistochemical analysis of NEK9 (Ph...



Immunofluorescent analysis of NEK9 (Phos...