



Product Datasheet

PAK1/2/3 antibody (orb315643)


Description	Rabbit polyclonal antibody to PAK1/2/3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human PAK1/2/3. The exact sequence is proprietary.
Target	PAK1; PAK2; PAK3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q61036 , Q13153 , Q62829 , O75914 , Q64303 , Q8CIN4 , O88643 , P35465 , Q13177
Entrez	5063 , 29431 , 100910732 , 224105 , 18481 , 29433 , 5062 , 5058
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of PAK1/2/3 expres...



Immunohistochemical analysis of PAK1/2/3...



Immunofluorescent analysis of PAK1/2/3 S...