



Product Datasheet

MYPT1 (Phospho-T853) antibody (orb315657)

Description	Rabbit polyclonal antibody to MYPT1 (Phospho-T853)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T853 of human MYPT1 protein. The exact sequence is proprietary.
Target	PPP1R12A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q14974, Q9DBR7, Q10728
Entrez	116670, 4659, 17931
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MYPT1 (Phospho-...



Immunohistochemical analysis of MYPT1 (P...