

Product Datasheet

PKA C alpha/beta/gamma (Phospho-T198) antibody (orb315659)

Description

Rabbit polyclonal antibody to PKA C

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat, Sheep

Conjugation

Unconjugated

Tested Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T198 of human PKA C alpha/beta/gamma protein. The exact sequence is proprietary.

Target

PRKACA; PRKACB

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P68182](#), [P27791](#), [P22694](#), [P68181](#), [P05132](#), [P22612](#), [P17612](#)

Entrez

[293508](#), [5567](#), [5566](#), [18747](#), [18749](#)


Dilution Range

WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500


Expiration Date

12 months from date of receipt.


kDa A B C D E F G



Western blot analysis of PKA C alpha/beta/gamma protein.



Immunohistochemical analysis of PKA C alpha/beta/gamma protein.



Immunofluorescent analysis of PKA C alpha/beta/gamma protein.