



Product Datasheet

DNA-PKcs (Phospho-T2647) antibody (orb315668)


Description	Rabbit polyclonal antibody to DNA-PKcs (Phospho-T2647)
Species/Host	Rabbit
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T2647 of human DNA-PKcs protein. The exact sequence is proprietary.
Target	PRKDC
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P78527
Entrez	5591
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of DNA-PKcs (Phospho-T2647) antibody (orb315668).



Immunohistochemical analysis of DNA-PKcs (Phospho-T2647) antibody (orb315668).



Immunofluorescent analysis of DNA-PKcs (Phospho-T2647) antibody (orb315668).