



Product Datasheet

RPA2 (Phospho-T21) antibody (orb315689)


Description	Rabbit polyclonal antibody to RPA2 (Phospho-T21)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T21 of human RPA2 protein. The exact sequence is proprietary.
Target	RPA2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q63528, P15927
Entrez	59102, 6118
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of RPA2 (Phospho-T21)



Immunohistochemical analysis of RPA2 (Phospho-T21)



Immunofluorescent analysis of RPA2 (Phospho-T21)