



Product Datasheet

Topoisomerase 2 alpha antibody (orb315729)

Description	Rabbit polyclonal antibody to Topoisomerase 2 alpha
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Topoisomerase 2 alpha. The exact sequence is proprietary.
Target	TOP2A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q01320 , P11388
Entrez	7153 , 21973
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Topoisomerase 2...



Immunohistochemical analysis of Topoisom...