



Product Datasheet

p53 (Phospho-S20) antibody (orb315733)


Description	Rabbit polyclonal antibody to p53 (Phospho-S20)
Species/Host	Rabbit
Reactivity	Human, Porcine, Primate, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S20 of human p53 protein. The exact sequence is proprietary.
Target	TP53
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P04637
Entrez	7157
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of p53 (Phospho-S20) antibody.



Immunohistochemical analysis of p53 (Phospho-S20) antibody.



Immunofluorescent analysis of p53 (Phospho-S20) antibody.