



## Product Datasheet

### TRAPPC3 antibody (orb315735)

|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TRAPPC3                                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Tested Applications</b> | IH, WB                                                                                                                                 |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TRAPPC3. The exact sequence is proprietary. |
| <b>Target</b>              | TRAPPC3                                                                                                                                |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                              |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                              |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                       |
| <b>Note</b>                | For research use only                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Source</b>              | Rabbit                                                                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q5U1Z2</a> , <a href="#">O55013</a> , <a href="#">O43617</a>                                                               |
| <b>Entrez</b>              | <a href="#">362599</a> , <a href="#">27096</a> , <a href="#">27095</a>                                                                 |
| <b>Dilution Range</b>      | WB: 1-500:1000, IHC-P: 1-100:200                                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                        |



Western blot analysis of TRAPPC3 expression...



Immunohistochemical analysis of TRAPPC3

...