



Product Datasheet

VAMP1/2/3 antibody (orb315739)

Description	Rabbit polyclonal antibody to VAMP1/2/3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human VAMP1/2/3. The exact sequence is proprietary.
Target	VAMP1; VAMP2; VAMP3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P63025 , P23763 , P63024 , P63045 , Q62442 , P63027 , Q15836 , P63044 , Q63666
Entrez	25624 , 9341 , 22317 , 6843 , 29528 , 6844 , 24803 , 22319 , 22318
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of VAMP1/2/3 expression.



Immunohistochemical analysis of VAMP1/2/3 expression.