



Product Datasheet

TAK1 (Phospho-T187) antibody (orb315801)


Description	Rabbit polyclonal antibody to TAK1 (Phospho-T187)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T187 of human TAK1 protein. The exact sequence is proprietary.
Target	MAP3K7
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q62073, O43318, P0C8E4
Entrez	100910771, 6885, 26409, 313121
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of TAK1 (Phospho-T187)



Immunohistochemical analysis of TAK1 (Phospho-T187)



Direct ELISA antibody dose-response curve