



Product Datasheet

NT5C1A antibody (orb315806)


Description	Rabbit polyclonal antibody to NT5C1A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human NT5C1A. The exact sequence is proprietary.
Target	NT5C1A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	A3KFX0 , Q9BXI3
Entrez	230718 , 84618
Dilution Range	WB: 1-500:1000, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of NT5C1A expressi...



Immunohistochemical analysis of NT5C1A S...



Immunofluorescent analysis of NT5C1A sta...