


Product Datasheet

TAL1/2 (AcK221+222/K36+37) antibody (orb315818)

Description	Rabbit polyclonal antibody to TAL1/2 (Acetyl-K221+222/K36+37)
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K221+222/K36+37 of human TAL1/2 protein. The exact sequence is proprietary.
Target	TAL1; TAL2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q16559 , P22091 , P17542 , Q62282
Entrez	21350 , 6887 , 21349 , 6886
Dilution Range	WB: 1-500:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of TAL1/2 (Acetyl-...