


Product Datasheet


14-3-3 theta/tau (Phospho-S232) antibody (orb315829)

Description


Rabbit polyclonal antibody to 14-3-3



Western blot analysis of 14-3-3 theta/ta...



Immunohistochemical analysis of 14-3-3 t...



Immunofluorescent analysis of 14-3-3 the...

Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S232 of human 14-3-3 theta/tau protein. The exact sequence is proprietary.
Target	YWHAQ
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P27348 , P68255 , P68254
Entrez	10971 , 25577 , 22630 , 102634437
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.