

Product Datasheet

TAK1 (Phospho-T184) peptide (orb315863)

Description	The peptide is used to TAK1 (pT184) antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAP3K7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	POC8E4, Q62073, O43318
Entrez	313121, 6885, 100910771, 26409
Expiration Date	12 months from date of receipt.