

Product Datasheet

ERK1/2 peptide (orb315864)

Description	ERK1/2 peptide
Reactivity	Human, Mouse, Zebrafish
Tested Applications	Blocking
Target	MAPK1; MAPK3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63844 , P63085 , P27361 , P28482
Entrez	26413 , 26417 , 5595 , 5594
Expiration Date	12 months from date of receipt.