

Product Datasheet

PDCD4 (Phospho-S67) peptide (orb315912)

Description	The peptide is used to PDCD4 (pS67) antibody
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PDCD4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61823 , Q53EL6 , Q9JID1
Entrez	64031 , 27250 , 18569
Expiration Date	12 months from date of receipt.