

Product Datasheet

SNAI1/2 (Phospho-S246) peptide (orb315974)

Description	The peptide is used to SNAI1/2 (pS246) antibody
Reactivity	Human, Mouse, Primate
Tested Applications	Blocking
Target	SNAI1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O95863, Q02085
Entrez	6615, 20613
Expiration Date	12 months from date of receipt.