

Product Datasheet

FGF13 peptide (orb316047)

Description	The peptide is used to FGF13 antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FGF13
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70377, Q92913, Q9ERW3
Entrez	14168, 2258, 84488
Expiration Date	12 months from date of receipt.