

Product Datasheet

ERK1/2 (Phospho-T202/Y204) peptide (orb316068)

Description	ERK1/2 (Phospho-T202/Y204) peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAPK1; MAPK3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P21708, P63085, P63086, Q63844, P27361, P28482
Entrez	50689, 116590, 26413, 5594, 26417, 5595
Expiration Date	12 months from date of receipt.