

Product Datasheet

CTAG1A antibody (orb31609)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CTAG1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17992
Uniprot ID	P78358
NCBI	NP_640343.1, NP_001318.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

