

Product Datasheet

XRN2 peptide (orb316094)

Description	The peptide is used to XRN2 antibody
Reactivity	Human, Mouse, Primate, Zebrafish
Tested Applications	Blocking
Target	XRN2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H0D6, Q9DBR1
Entrez	24128, 22803
Expiration Date	12 months from date of receipt.