


Product Datasheet

Human SHBG / Sex Hormone Binding Globulin ELISA Kit (orb316527)

Description	Human SHBG / Sex Hormone Binding Globulin ELISA Kit
Reactivity	Human
Range	0.031 nmol/L - 2 nmol/L
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L30-H402
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	43779 MW
Uniprot ID	P04278
Sensitivity	0.001 nmol/L
Dilution Range	For serum, dilution ratio is 1:1000 and concentration is around 130nmol/l. For plasma, dilution ratio is 1:10 and concentration is around 2nmol/l.
Expiration Date	Please enquire.



Human
SHBG
PicoKine
ELISA Kit
standard
C...