


Product Datasheet

Rat CD44 ELISA Kit (orb316533)

Description	Rat CD44 ELISA Kit
Reactivity	Rat
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q22-E271
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	55946 MW
Uniprot ID	P26051
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 10-25 ng/ml.
Expiration Date	Please enquire.



Rat CD44
PicoKine
ELISA Kit
standard
cur...