

Product Datasheet

Human CD73/NT5E/Ecto-5'-nucleotidase ELISA Kit (orb316537)

Description Human CD73/NT5E/Ecto-5'-nucleotidase ELISA Kit

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: CHO; Immunogen sequence: W27-K547

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).


MW 63368 MW

Uniprot ID [Q6NZX3](#)

Sensitivity 10 pg/ml

Dilution Range Dilution ratio of 1:10, concentration in serum and plasma is around 6.4 ng/ml.

Expiration Date Please enquire.



Human
CD73
PicoKine
ELISA Kit
standard
C...