

Product Datasheet

GRK3/ADRBK2 Antibody (orb316545)

Description

GRK3/ADRBK2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IF, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human GRK3 (635-669aa ESDPEFVQWKELNETFKEAQRLLRRAPKFLNKPR), different from the related mouse and rat sequences by five amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

79710 MW

Uniprot ID


P35626

Dilution Range

Western blot, 0.1-0.5µg/ml, Human
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human
Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Flow
Cytometry
analysis of
A431 cells
us...



WB analysis of
GRK using
anti-GRK
antibo...



IF analysis of
GRK3 using
anti-GRK3
anti...