

Product Datasheet

B3GNT8 Antibody (orb316549)

Description

B3GNT8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human B3GNT8 (360-397aa ADRTADHCAFRNLLLVRPLGPQASIRLWKQLQDPRLQC), different from the related mouse sequence by sixteen amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

43396 MW

Uniprot ID

[Q7Z7M8](#)

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.

130KD —

100KD —


70KD —

55KD —

35KD —

25KD —

WB analysis of B3GNT8 using anti-B3GNT8 ...



IHC analysis of B3GNT8 using anti-B3GNT8...