

Product Datasheet

BMP6 Antibody (orb316553)

Description

BMP6 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human BMP6 (383-419aa QSRNRSTQSQDVARVSSASDYNSELKTACRKHELYV), different from the related mouse sequence by six amino acids, and from the related rat sequence by three amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

57226 MW

Uniprot ID


P22004

Dilution Range

Western blot, 0.1-0.5µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human ELISA , 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



Flow
Cytometry
analysis of
HeLa cells us...



WB analysis of
BMP6 using
anti-BMP6
anti...