


Product Datasheet

MIG/CXCL9 Antibody (orb316557)

Description	MIG/CXCL9 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	E. coli-derived human CXCL9 recombinant protein (Position: T23-T125). Human CXCL9 shares 68.2% amino acid (aa) sequence identity with mouse CXCL9.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14019 MW
Uniprot ID	Q07325
Dilution Range	ELISA , 0.1-0.5µg/ml, Human, -Western blot, 0.1-0.5µg/ml, Human
Expiration Date	12 months from date of receipt.



WB analysis of CXCL9 using anti-CXCL9 an...