

Product Datasheet

FAM102B antibody (FITC) (orb31678)

Description	FAM102B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein FAM102B protein (92-360AA)
Target	FAM102B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T8I3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.