


Product Datasheet

HAO2 antibody (orb316839)

Description	Goat polyclonal antibody to HAO2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HAO2 (aa54-67)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 38kDa band observed in Human Kidney lysates (calculated MW of 38.8kDa according to NP_057611.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	DVSEVDTRTTIQGE
Clonality	Polyclonal
MW	38.8
Entrez	51179
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...