

Product Datasheet

LIMS2 antibody (orb31684)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target LIMS2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 38916

Uniprot ID [Q7Z4I7](#)

NCBI [NP_001129509.2](#), [NP_001154876.1](#),
[NP_001243471.1](#), [NP_060450.2](#), [NP_001154875.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

