

Product Datasheet

FNDC7 antibody (orb31691)

Description	FNDC7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Fibronectin type III domain-containing protein 7 protein (234-533AA)
Target	FNDC7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5VTL7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.