

Product Datasheet

AMACR / p504S Antibody (Prostate Cancer Marker) (orb316989)

Description

This antibody recognizes a protein of 44kDa,



p504S. It i...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Other

Tested

IHC-P, WB

Applications

Immunogen

Full length human recombinant protein was used as the immunogen for the p504S antibody.

Preservatives

Bioreactor concentrate with 0.05% sodium azide

Storage

Store the p504S antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the p504S antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 10-20 min followed by cooling at RT for 20 min2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

Bioreactor concentrate with 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Monoclonal

Purity

Unpurified high titer supernatant

Uniprot ID

Q9UHK6

Hazard

Information


This p504S antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1:50-1:100 for 30-60 min at RT, Western blot: 1:-50-1:100

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human pr...



Immunohistochemical staining of Human pr...