

Product Datasheet

ATG1 antibody (orb31721)

Description	ATG1 antibody
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Candida glabrata Serine/threonine-protein kinase ATG1 protein (11-312AA)
Target	ATG1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6FL58
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.