

Product Datasheet

CDw75 Antibody (orb317313)

Description

Recognizes a neuraminidase-sensitive sialoprotein

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, IF, IHC-P

Applications

Immunogen

Nuclei from pokeweed mitogen-stimulated peripheral blood lymphocytes were used as the immunogen for the CDw75 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CDw75 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the CDw75 antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgM, kappa

Clonality

Monoclonal

Purity

PEG precipitation

Uniprot ID

P15907

Hazard

Information


This CDw75 antibody is available for research use only.

Dilution Range


Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml, Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human sp...



Flow cytometric analysis of human lympho...