

Product Datasheet

ATG8 antibody (FITC) (orb31751)

Description	ATG8 antibody (FITC)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	FITC
Immunogen	Recombinant Candida glabrata Autophagy-related protein 8 protein (1-116AA)
Target	ATG8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6FXR8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.