

Product Datasheet

APOC2 antibody (orb317549)

Description	APOC2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	KLH conjugated synthetic peptide derived from human APOC2 (51-97/97 aa)
Target	APOC2
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9 kDa
Uniprot ID	P02655
Expiration Date	12 months from date of receipt.