


Product Datasheet

ADCY3 antibody (orb317553)

Description	ADCY3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human ADCY3 (651-750/1144 aa)
Target	ADCY3
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	126 kDa
Uniprot ID	O60266
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of
Lane 1:
(mouse/spleen)