

Product Datasheet

ARP3 (phospho-Ser418) antibody (orb317559)

Description	ARP3 (phospho-Ser418) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human ARP3 around the phosphorylation site of Ser418 (VM(p-S))
Target	phospho-ARP3 (Ser418)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	P61158
Expiration Date	12 months from date of receipt.