


Product Datasheet

CXCR4 antibody (orb317585)

Description	CXCR4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse CXCR4 (151-300/359 aa)
Target	CXCR4
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P61073
Dilution Range	WB=1:500-2000, Flow-Cyt=3µg/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Mouse
spleen...