

Product Datasheet

acpS antibody (HRP) (orb31759)

Description	acpS antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus Holo-[acyl-carrier-protein] synthase protein (1-119AA)
Target	acpS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6GF02
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.