

Product Datasheet

CDC37 (phospho-Ser13) antibody (orb317591)

Description	CDC37 (phospho-Ser13) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CDC37 around the phosphorylation site of Ser13 (EV(p-S)DD)
Target	phospho-CDC37 (Ser13)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q16543
Expiration Date	12 months from date of receipt.