


Product Datasheet

CD68 antibody (orb317598)

Description	CD68 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse CD68 (151-250/335 aa)
Target	CD68
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P34810
Dilution Range	WB=1:500-2000, IF=1:100-500, Flow-Cyt=1µg/Test, IHC-P=1:100-500, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Rsc96 cell u...