


Product Datasheet

CD68 antibody (orb317599)

Description	CD68 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from mouse CD68 (111-210/335 aa)
Target	CD68
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P34810
Dilution Range	Flow-Cyt=3ug/Test
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Brain Lys...