


Product Datasheet

EAAT2 antibody (orb317613)

Description	EAAT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human EAAT2 (351-450/574 aa)
Target	EAAT2
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62 kDa
Uniprot ID	P43006
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ELISA=1:5000-10000, IHC-F=1:100-500, ICC=1:100, IF=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brai...