



Product Datasheet

MYPT1 (Phospho-Thr696) antibody (orb317659)

Description	MYPT1 (Phospho-Thr696) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MYPT1 around the phosphorylation site of Thr696 (RS(p-T)QG)
Target	Phospho-MYPT1 (Thr696)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	O14974
Dilution Range	ICC=1:100, Flow-Cyt=1ug/Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of HepG2(black line).Prim...



ICC analysis of HepG2 cell using Phospho...