

Product Datasheet

RAB11FIP1 antibody (orb317772)

Description	RAB11FIP1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human RAB11FIP1 (1101-1200/1283 aa)
Target	RAB11FIP1
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	137 kDa
Uniprot ID	Q6WKZ4
Expiration Date	12 months from date of receipt.