



Product Datasheet

BRCA2 antibody (orb317801)


Description	BRCA2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Rad51 (81-180/339 aa)
Target	Rad51
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	Q06609
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:25, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human CO...



Western blot analysis of mouse Bone lysa...



Immunohistochemical staining of rat test...