


Product Datasheet

Integrin alpha-3 heavy chain antibody (orb317807)

Description	Integrin alpha-3 heavy chain antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Integrin alpha-3 heavy chain (33-100/1051 aa)
Target	Integrin alpha-3 heavy chain
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93/114 kDa
Uniprot ID	P26006
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Brain tis...